




Hospira Medical UK Ltd

Company Profile

Hospira Medical UK Ltd. Specializes in developing and marketing of Covid-19 Medical Protective Equipment, PPE, oxygen concentrator and other healthcare related products. All of our products are FDA/CE certified in accordance with US/EU certification. We also have a perfect quality management system, excellent sales services, and offer our clients pleasant shipping services. Hospira UK is a subsidiary of Pfizer, the leading manufacturer of Pfizer-BioNTech COVID-19 vaccine.

CONTENTS

Part 1 / Medical Mask & Respirators



Surgical Face Mask Features



【Standard】 YY0469-2011, EN14683 Type IIR

【Certificate】 CE Certificate directive 93/42/EEC

Have been met for the products, Issued by Notified body:
TUV Rheinland

【Sterile】 Yes

【Specification】 17*9.5cm, 28gsm pp+25gsm Melt-
blown+25gsm PP

【Package】 50pcs/box, 500pcs/carton

【Package Size】 : 17.5*17cm

【Carton Size】 : 52*39*52cm

【Particulate Filtration Efficiency】 :

>98% **【Splash Resistance**

Pressure】 :120mmHg **【Color】** : Blue

【Validity】 : 2 years



Model Number: ZH3161

Name: Folding Mask **Class:** FFP2 NR

Valve: No

Color: White / Blue

Layer: 4 ply

Material: 29% PP Spunbond, 40% PP Melt blown
Fabric, 31% Polyester fiber

Packaging: 30PCS/BOX 20BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3161V FFP2NR

Name: Folding Mask **Class:** FFP2 NR

Valve: No

Color: White / Blue

Layer: 4 ply

Material: 29% PP Spunbond, 40% PP Melt blown
Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3310 FFP2 NR

Name: Folding Cup Shape Mask with Valve

Class: FFP2 NR

Valve: Yes

Color: White / Red

Layer: 3 ply

Material: 29% PP Spunbond, 40% PP Melt
blown Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3310V FFP2 NR

Name: Folding Cup Shape Mask with Valve

Class: FFP2 NR

Valve: Yes

Color: White / Red

Layer: 3 ply

Material: 29% PP Spunbond, 40% PP Melt blown

Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3361 FFP3NR
Name: Folding Mask **Class:** FFP3 NR

Valve: No

Color: White / Blue

Layer: 4 ply

Material: 29% PP Spunbond, 40% PP Melt blown
Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3361V FFP3NR

Name: Folding Mask **Class:** FFP3 NR

Valve: No

Color: White / Blue

Layer: 4 ply

Material: 29% PP Spunbond, 40% PP Melt
blown Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3081 FFP3 NR

Name: Folding Cup Shape Mask with Valve

Class: FFP3 NR

Valve: Yes

Color: White / Red

Layer: 3 ply

Material: 29% PP Spunbond, 40% PP Melt blown Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg



Model Number: ZH3081V FFP3NR

Name: Folding Cup Shape Mask with Valve

Class: FFP3 NR

Valve: Yes

Color: White / Red

Layer: 3 ply

Material: 29% PP Spunbond, 40% PP Melt blown
Fabric, 31% Polyester fiber

Packaging: 20PCS/BOX 10BOX/CTN

Dimension: 69*28.5*34cm

G.W: 6.5kg

Hospira Medical UK Ltd



Website: www.hospiramedical.co.uk

E-mail: info@hospiramedical.co.uk

Mobile/Whatsapp: +44 192 695 5390